Introduction to ImageJ

Improve your image analysis pipeline

2017-04-24

Xavier Draye and Guillaume Lobet

Image analysis pipeline

Image analysis is only one of the step in the global data analysis pipeline:


- 1. Image acquisition
- 2. Image analysis
- 3. Data analysis

Sometimes it is easier to improve previous or future steps rather than spend too much time on the image analysis in itself...

Improve image acquisition

- Better lighting conditions
- Keep same disctance between camera and object
- Increase contrast between background and object
- Reduce the number of artefact in the object
- Try different sensors (visible, infrared, multispectral)

Improve the image acquisition




A. Camera
B. Fixed distance
C. Black background
D. Scale

Improve image analysis

- Calibrate your pipeline
- Find a robust ground-truth dataset
 - manual measurements
 - simulated images
- Estimate the error

Estimate the error



Improve data processing

- Filter out noisy data
- Correctly identify outlyers